

## Kompatibilita manžiet OMRON

Obchodné označenie tlakomeru	Modelové označenie tlakomeru	Na obvod paže (v cm)					
		17 - 22	22 - 32	22 - 42			
		mäkká manžeta prevlečená cez sponu				pevne tvarovaná manžeta	
		OMRON <u>CS</u>	OMRON <u>CM</u>	OMRON <u>CW</u>	OMRON <u>Easy</u>	OMRON <u>Comfort</u>	OMRON <u>Intelli</u>
<b><u>M2 Basic</u></b>	<b><u>HEM-7121J</u></b>	○	●	○	○		
M2 Basic	HEM-7120	○	●	○	○		
M2 Basic	HEM-7116	○	●	○			
M2 Compact	HEM-7102	○	●	○			
<b><u>M2</u></b>	<b><u>HEM-7121</u></b>	○	●	○	○		
M2	HEM-7119	○	●	○			
M2	HEM-7117	○	●	○			
<b><u>M3</u></b>	<b><u>HEM-7154</u></b>	○			●		
M3	HEM-7131	○	○		●		
M3 IT	HEM-7131U	○	○		●		
<b><u>M3 Comfort</u></b>	<b><u>HEM-7155</u></b>						●
M3 Comfort	HEM-7134						●
M3	HEM-7200	○	●	○			
M3W	HEM-7202			●			
M3 Intellisense	HEM-7051	○	●	○			
MX3 Plus	HEM-742	○	●	○			
<b><u>M4 Intelli IT</u></b>	<b><u>HEM-7155T</u></b>						●
M4-I	HEM-752		●	○			
M6 AC	HEM-7322	○			●		
M6W	HEM-7213	○		●			
M6	HEM-7211	○	●	○			
<b><u>M6 Comfort</u></b>	<b><u>HEM-7360</u></b>						●
M6 Comfort	HEM-7321						●
M6 Comfort	HEM-7223					●	
M6 Comfort	HEM-7221					●	
M6 Comfort IT	HEM-7322U						●
<b><u>M7 Intelli IT</u></b>	<b><u>HEM-7361</u></b>						●
M7 Intelli IT	HEM-7322T						●
M7	HEM-780					●	
M10-IT	HEM-7080IT					●	
i-C10	HEM-7070					●	
<b><u>M1 Compact</u></b>	<b><u>HEM-4022</u></b>	○	●	○			
M1	HEM-4030	○	●	○			
M1 Plus	HEM-4011C	○	●	○			
M1 Classic	HEM-442	○	●	○			

● manžeta v základnom balení tlakomeru

○ manžetu možno dokúpiť ako voliteľné príslušenstvo